


PATENT FIGURES AND ARCHITECTURE DIAGRAMS

FIG. 1: Overall System Architecture







FIG. 2: Hierarchical Status Computation Engine Components







FIG. 3: Cross-Entity Validation Framework Workflow



ToStatus: APPROVED
InitiatedBy: regulatory_reviewer@clinprecision.com




DEPENDENCY GRAPH BUILDER

Step 1: Identify Related Entities

PROTOCOL_VERSION: PV-456
├ STUDY: ST-123 (Status: UNDER_REVIEW)
├ SITES: [SI-789, SI-790, SI-791]
│ ├── SI-789 (Status: READY_FOR_ACTIVATION)
│ ├── SI-790 (Status: READY_FOR_ACTIVATION)
│ └ SI-791 (Status: UNDER_SETUP)
└ REGULATORY: REG-555 (Status: PENDING_APPROVAL)

Step 2: Build Real-time Dependency Graph



CONSTRAINT VALIDATOR

Business Rule Validation:

- ✓ Rule BR-001: Protocol Version can only be APPROVED if parent Study is not in TERMINATED status
Current: Study ST-123 = UNDER_REVIEW ✓
- ✓ Rule BR-002: Protocol Version approval requires at least one Site in READY_FOR_ACTIVATION status
Current: SI-789, SI-790 = READY_FOR_ACTIVATION ✓
- ⚠ Rule BR-003: All Sites should be ready before Protocol Version approval
Current: SI-791 = UNDER_SETUP ⚠ WARNING



REGULATORY COMPLIANCE CHECKER

FDA Guidelines Validation:





FIG. 4: Event-Driven State Synchronization Process



2. Event Classification

- └ Event Priority: HIGH (Protocol Version Approval)
- └ Impact Scope: MULTI_ENTITY
- └ Processing Mode: SYNCHRONOUS



IMPACT ANALYZER

Dependency Impact Analysis:

Primary Impacts:

STUDY: ST-123
Current Status: UNDER_REVIEW
Recommended Action: RECOMPUTE_STATUS
Reason: Active Protocol Version now available

Secondary Impacts:

SITES: [SI-789, SI-790, SI-791]
Current Status: READY_FOR_ACTIVATION
Recommended Action: ENABLE_ACTIVATION
Reason: Approved Protocol Version available

Tertiary Impacts:

PARTICIPANTS: All enrolled participants
Current Status: Various
Recommended Action: UPDATE_PROTOCOL_REFERENCE
Reason: Protocol Version reference update required



SYNCHRONIZATION ORCHESTRATOR

Transaction Management:

Phase 1: Preparation

- └ Create distributed transaction context
- └ Acquire locks on all affected entities
- └ Validate all proposed status changes

Phase 2: Execution

- └ Execute status updates in dependency order
- └ Update STUDY ST-123: UNDER_REVIEW → APPROVED
- └ Update SITES [SI-789, SI-790]: → READY_TO_ACTIVATE
- └ Update PARTICIPANTS: Protocol references





FIG. 5: Compliance Rule Engine Architecture



ICH-GCP 8.2: Quality assurance procedures



BUSINESS RULE REPOSITORY

Clinical Trial Workflow Rules:

- BR-001: Study cannot be ACTIVE without approved protocol
- BR-002: Sites require activated protocol for enrollment
- BR-003: Participants require active site for enrollment
- BR-004: Protocol versions follow sequential approval
- BR-005: Status changes require authorized user actions

Entity Relationship Rules:

- ER-001: One-to-many: Study to Protocol Versions
- ER-002: Many-to-many: Studies to Sites
- ER-003: Many-to-one: Participants to Sites
- ER-004: Many-to-one: Sites to Protocol Versions
- ER-005: Hierarchical: Study → Protocol → Site → Participant



RULE EXECUTION ENGINE

Rule Processing Pipeline:

1. Rule Selection
 - └ Context Analysis (Entity Type, Status Transition)
 - └ Applicable Rule Identification
 - └ Rule Priority Ordering
2. Rule Evaluation
 - └ Entity State Retrieval
 - └ Rule Condition Assessment
 - └ Compliance Status Determination
3. Rule Enforcement
 - └ Violation Prevention (Block invalid transitions)
 - └ Warning Generation (Flag potential issues)
 - └ Corrective Action Recommendation
4. Audit Trail Generation
 - └ Rule Execution Logging
 - └ Compliance Decision Recording
 - └ Regulatory Reporting Preparation

COMPLIANCE ASSESSMENT RESULT

ComplianceStatus: COMPLIANT
RulesEvaluated: 12
RulesPassed: 11
RulesFailed: 0
WarningsGenerated: 1
WarningDetails: "Site readiness optimization recommended"
RegulatoryReporting: AUTO_GENERATED
AuditTrailId: AT-COMP-789123
NextReviewDate: 2025-10-28