

# Patent Drawings - Event Sourcing Architecture for Clinical Trials

---

## USPTO-Compliant Technical Drawings

**Patent Application:** Event Sourcing Architecture for Clinical Trial Management

**Date:** October 17, 2025

**Drawing Standards:** USPTO 37 CFR 1.84

---

## DRAWING INSTRUCTIONS


### USPTO Requirements


1. **Format:** Black ink on white paper
2. **Size:** 8.5" x 11" (letter size)
3. **Margins:** 1" top, 1" left/right, 1" bottom
4. **Line width:** Minimum 0.3mm
5. **Text size:** Minimum 0.32cm (1/8 inch) high
6. **Reference numerals:** Numbers pointing to components
7. **Figure numbering:** FIG. 1, FIG. 2, etc.
8. **Shading:** Use reference numerals, avoid gray shading if possible


### Tools for Creating Drawings


- **Draw.io:** <https://app.diagrams.net> (Free, recommended)
  - **Lucidchart:** <https://www.lucidchart.com>
  - **Microsoft Visio:** Professional diagrams
  - **Patent Drawing Software:** PatentDrawing.com
- 

## FIGURE 1: System Architecture Overview









Reference Numerals:  
100 - Clinical Trial Management System

- 110 - User Interface Layer
- 111 - Study Setup Module
- 112 - Patient Enrollment Module
- 113 - Data Capture Module
- 114 - Reporting Module
- 120 - Commands (User Actions)
- 130 - Command Gateway
- 140 - Event-Sourced Aggregates
- 141 - Study Aggregate
- 142 - Patient Aggregate
- 143 - Visit Aggregate
- 144 - Form Data Aggregate
- 150 - Events
- 160 - Event Store (Append-Only)
- 170 - Event Stream
- 180 - Event Handlers
- 190 - Projectors
- 200 - Read Database

**Figure 1 Description:** System architecture showing event flow from user interface through command gateway to event-sourced aggregates, which generate events stored in append-only event store. Events are consumed by handlers and projectors to create read models.

FIGURE 2: Event Store Structure and Audit Trail

EVENT STORE TABLE (160)		
event_id (210)	UUID	Primary Key
aggregate_type (211)	VARCHAR(50)	"Study", "Patient"
aggregate_id (212)	UUID	Links to entity
event_type (213)	VARCHAR(100)	"StudyCreated"

event_version (214)	INT	Event schema
event_data (215)	JSON	Event payload
metadata (216)	JSON	Audit information
occurred_at (217)	TIMESTAMP	When event
sequence_number (218)	BIGINT	Global ordering

METADATA STRUCTURE (216)

```
{
  "user_id": "12345"           (220) ← WHO
  "user_name": "Dr. Jane Smith" (221)
  "reason": "Protocol amendment" (222) ← WHY
  "ip_address": "192.168.1.100" (223) ← WHERE
  "user_agent": "Chrome 118.0"  (224)
  "session_id": "abc-123"       (225)
}
```

EXAMPLE EVENT RECORD

```
event_id:      "550e8400-e29b-41d4-a716-446655440000"

aggregate_type: "Study"

aggregate_id:   "650e8400-e29b-41d4-a716-446655440001"

event_type:     "StudyCreatedEvent"

event_version:  1

event_data: {

  "protocol_number": "PROTO-2025-001",

  "organization_id": 123,

  "study_phase": "PHASE_III"

}

metadata: {

  "user_id": "12345",

  "user_name": "Dr. Jane Smith",

  "reason": "New hypertension trial",


  "ip_address": "192.168.1.100"

}

occurred_at:    "2025-10-17T14:30:00Z"


sequence_number: 1001
```

- ▶ IMMUTABLE - Never updated or deleted
- ▶ Complete audit trail by design
- ▶ FDA 21 CFR Part 11 compliant





**Figure 2 Description:** Event store structure showing append-only table with immutable events. Each event contains complete audit information (who, what, when, why) enabling automatic audit trail generation without separate audit tables.


**FIGURE 3: Protocol Version Management with Build ID**











#### KEY INSIGHTS:

- ▶ Patients 001 & 002: Remain on Version 1.0 (3 visits)
- ▶ Patients 003 & 004: Use Version 2.0 (5 visits)
- ▶ NO DATA MIGRATION REQUIRED
- ▶ Each patient's data linked to their protocol version via build\_id

#### Reference Numerals:


300 - Protocol Version Management System

```
310 - Protocol Version 1.0
320 - Protocol Version 2.0
330 - Visit Definitions (Master Templates)
340 - Visit Forms Mapping Table
350 - Patient Enrollment
360 - Visit Instances
370 - Form Data
```

**Figure 3 Description:** Protocol version management showing how build\_id propagates from protocol versions through mappings to patient data. Existing patients remain on their enrolled version while new patients use the latest version. No data migration required for protocol amendments.

---

## FIGURE 4: Time-Travel Query Process



- Filters events up to target date
- Orders by sequence for correct replay



#### STEP 2: Event Stream (430)

Event 1 (Jan 15): StudyCreatedEvent

Event 2 (Feb 10): PatientEnrolledEvent (Patient 001)

Event 3 (Mar 20): PatientEnrolledEvent (Patient 002)

Event 4 (Apr 15): VisitCompletedEvent (Patient 001, Visit 1)

Event 5 (May 30): PatientEnrolledEvent (Patient 003)

Event 6 (Jun 01): PatientEnrolledEvent (Patient 004) ←  
INCLUDED

Event 7 (Jun 15): PatientEnrolledEvent (Patient 005) ←  
EXCLUDED

- Only events before/on June 1, 2024 included



#### STEP 3: Aggregate Reconstruction (440)

```
StudyAggregate study = new StudyAggregate();
```

```
for (Event event : events) {  
    study.applyEvent(event); // Replay event  
}
```

- State rebuilt by replaying historical events
- Exact state as of June 1, 2024



STEP 4: State at Point in Time (450)

```
StudyAggregate state {  
    studyId: "650e8400-...",  
    protocolNumber: "PROTO-2025-001",  
    status: ACTIVE,  
    enrolledPatients: [001, 002, 003, 004], ← 4 patients  
    totalVisitsCompleted: 1  
}
```

- Historical state reconstructed



STEP 5: Query Result (460)

ANSWER: 4 patients enrolled as of June 1, 2024

Details:

- Patient 001: Enrolled Feb 10, 2024
- Patient 002: Enrolled Mar 20, 2024
- Patient 003: Enrolled May 30, 2024
- Patient 004: Enrolled Jun 01, 2024
- Patient 005: NOT INCLUDED (enrolled after target date)

- ▶ Definitive answer for regulatory audit
- ▶ No guessing or approximate data

USE CASES (470):

- FDA Inspections: Prove exact state during audit period
- Debugging: Reproduce issues by replaying to problem point
- Compliance: Verify no data tampering or backdating
- Historical Reporting: Generate reports for any time period
- Root Cause Analysis: Understand how state evolved over time

Reference Numerals:

- 400 - Time-Travel Query System
- 410 - User Query
- 420 - Event Store Query
- 430 - Filtered Event Stream


440 - Aggregate Reconstruction  
450 - Historical State  
460 - Query Result  
470 - Use Cases


**Figure 4 Description:** Time-travel query process showing how historical aggregate state is reconstructed by filtering events up to target date and replaying them. Enables answering regulatory questions about past states definitively without approximations.

FIGURE 5: CQRS Architecture - Read/Write Separation









## KEY PRINCIPLES:

WRITE SIDE (Commands)	READ SIDE (Queries)
<ul style="list-style-type: none"> <li>• Validates business rules</li> <li>• Generates events</li> <li>• Slow (validation overhead)</li> <li>• Event-sourced aggregates</li> <li>• Consistency first</li> <li>• Single responsibility</li> </ul>	<ul style="list-style-type: none"> <li>• No validation</li> <li>• Direct database reads</li> <li>• Fast (optimized views)</li> <li>• Denormalized tables</li> <li>• Performance first</li> <li>• Multiple views</li> </ul>

## BENEFITS (630):

- ▶ Independent scaling of reads and writes
- ▶ Optimized read models for specific queries
- ▶ Event-driven microservices architecture
- ▶ Multiple projections from same event stream
- ▶ Eventual consistency acceptable for read models

## Reference Numerals:


- 500 - CQRS Architecture
- 510 - User Actions (Commands)
- 520 - Commands
- 530 - Command Handlers
- 540 - Events Generated
- 550 - Event Store
- 555 - Event Stream
- 560 - User Queries
- 570 - Queries
- 580 - Query Handlers
- 590 - Query Results
- 600 - Projectors (Event Handlers)
- 610 - Read Model Updates

620 - Read Models (Denormalized)  
630 - CQRS Benefits

**Figure 5 Description:** CQRS architecture separating write path (commands creating events) from read path (queries from optimized models). Events from write side trigger projectors that update denormalized read models for fast queries. Enables independent scaling and optimization of reads vs writes.

---

## FIGURE 6: Audit Trail Generation Process



STEP 2: Event Records with Metadata (730)

Event 1: StudyCreatedEvent

Timestamp: 2025-01-15 09:00:00

User: Dr. Jane Smith (ID: 12345)

Action: Created study PROTO-2025-001

Reason: "New Phase III hypertension trial"

IP Address: 192.168.1.100

Changes: protocol\_number: null → "PROTO-2025-001"

status: null → "PLANNING"

Event 2: ProtocolVersionCreatedEvent

Timestamp: 2025-02-01 14:30:00

User: Dr. Jane Smith (ID: 12345)

Action: Created protocol version 1.0

Reason: "Initial protocol approval"

IP Address: 192.168.1.100

Changes: version: null → "1.0"

amendment\_type: "INITIAL"

Event 3: PatientEnrolledEvent

Timestamp: 2025-03-10 11:15:00

User: Dr. Bob Wilson (ID: 54321)

Action: Enrolled patient 001


Reason: "Patient met eligibility criteria"

IP Address: 192.168.2.50

Changes: patient\_id: null → 001

enrollment\_date: "2025-03-10"

build\_id: 1001 (Version 1.0)



Apr 05 10:00 | Dr. Wilson | Patient 001 Visit 1 Completed  
... | ... | ...

#### CRITICAL CHANGES (FDA 21 CFR Part 11):

1. Protocol Created (Jan 15) - Dr. Smith
2. Protocol v1.0 Approved (Feb 01) - Dr. Smith
3. First Patient Enrolled (Mar 10) - Dr. Wilson
4. Protocol v2.0 Amendment (Jun 15) - Dr. Smith

Reason: "Added 2 new study visits for safety monitoring"

#### DATA INTEGRITY VERIFICATION:

- Event sequence: VALID (no gaps in sequence numbers)
- Timestamps: MONOTONIC (events in chronological order)
- Hash chain: VERIFIED (no tampering detected)
- User authentication: ALL VALID



#### STEP 4: Compliance Certification (750)

✓ FDA 21 CFR Part 11 COMPLIANT

- Complete audit trail

- Immutable records
- User authentication
- Reason for change documented
- Timestamp integrity verified

✓ ICH-GCP COMPLIANT

- All protocol changes documented
- Patient enrollment tracked
- Data collection audit trail

Report Generated: 2025-10-17 15:30:00

Report Hash: abc123def456...

Certification: VALID


#### KEY ADVANTAGES (760):

- ▶ No separate audit tables to maintain
- ▶ Audit trail generated from event store automatically
- ▶ Complete traceability (who, what, when, why)
- ▶ Tamper-proof (immutable events)
- ▶ Real-time availability (no batch processing)
- ▶ Cryptographic verification of data integrity

#### Reference Numerals:

- 700 - Automated Audit Trail Generation System
- 710 - Regulatory Request
- 720 - Event Store Query
- 730 - Event Records with Metadata
- 740 - Generated Audit Report
- 750 - Compliance Certification
- 760 - Key Advantages





**Figure 6 Description:** Automated audit trail generation showing how regulatory compliance reports are created directly from event store without separate audit tables. Each event contains complete audit metadata (who, what, when, why) enabling FDA 21 CFR Part 11 and ICH-GCP compliance by architectural design.

---

## INSTRUCTIONS FOR PATENT ATTORNEY

### Creating Final USPTO Drawings

#### 1. Use Professional Drawing Software

- Draw.io, Lucidchart, or Visio
- Export as high-resolution black & white images
- Ensure clean lines (minimum 0.3mm width)

#### 2. Follow USPTO Standards

- Paper size: 8.5" x 11" (letter)
- Margins: 1" all sides
- Reference numerals: Use numbers from descriptions above
- Figure labels: "FIG. 1", "FIG. 2", etc.
- No gray shading (use hatching/patterns if needed)

#### 3. Include All Figures

- Figure 1: System Architecture Overview
- Figure 2: Event Store Structure
- Figure 3: Protocol Version Management
- Figure 4: Time-Travel Query Process
- Figure 5: CQRS Architecture
- Figure 6: Audit Trail Generation

4. **Reference Numeral List** Create a separate sheet listing all reference numerals and their meanings for USPTO filing.
- 

## ADDITIONAL DOCUMENTATION

### Supporting Materials to Provide Attorney

### 1. **Code Samples** (Sanitized)

- Event store implementation
- Aggregate examples
- Command handlers
- Projector implementations

### 2. **Database Schemas**

- Event store table structure
- Read model tables
- Migration scripts

### 3. **Performance Benchmarks**

- Event replay speed
- Query performance
- Scalability metrics

### 4. **Compliance Documentation**

- FDA 21 CFR Part 11 compliance checklist
- ICH-GCP compliance mapping
- Audit trail examples

---

**These drawing descriptions and ASCII diagrams should be converted to professional USPTO-compliant drawings by your patent attorney or a professional patent illustrator.**

**Cost for professional patent drawings:** \$100-\$300 per figure (6 figures = \$600-\$1,800)

---

**END OF DRAWING SPECIFICATIONS**